



VitroSure™ DNA FFPE Tissue Isolation Kit
(For 50 Preps)
 SKU#: VB-5001

Description:

The VitroSure™ DNA FFPE Tissue Isolation Kit delivers fast, reliable, and high-yield DNA extraction from formalin-fixed, paraffin-embedded (FFPE) tissue samples. Engineered for both genomic and mitochondrial DNA isolation, this kit features proprietary VitroSure™ DNA Elute Columns and advanced silica membrane technology for consistent purification results. With flexible elution volumes (20–100 µl) and optimized performance for small tissue samples, this cost-effective, ready-to-use kit is ideal for sensitive downstream applications including PCR, qPCR, NGS, and genotyping.

Technical Specifications:

Equipment needed	Microcentrifuge, heat block/bath (37°C, 56°C and 90°C)
DNA Type Isolated	Total DNA
Size Range	> 50 bp
Yield	Up to 25 µg total DNA can be eluted into ≥ 50 µl
Purity	Typical A260/A280 ≥ 1.8
Eluted DNA Storage	at ≤ -20°C
Sample Source	Tissue from paraffin block or tissue sections
Processing Capacity	FFPE Tissue : ≤25 mg or 2-8 sections at thickness of 7-10 µm with a surface area of 15-20 mm ²
Applicable For	PCR and Next generation sequencing (NGS), genotyping, Restriction enzyme digestion, SNP, etc.

Kit Contents:

Buffer VTL	15 ml
RNase A (100 mg/ml)	120 µl
Buffer VL	15 ml
Buffer VW1	30 ml
Buffer VW2	30 ml
Buffer VTE	6 ml
Proteinase K Powder	5 mg×4
Proteinase K Buffer	2 ml
VitroSure DNA Elute Columns	50
Collection Tubes (2ml)	100

Storage

Store the Proteinase K Powder and RNase A (100 mg/ml) at -20°C. After reconstitution of Proteinase K, store the solution at -20°C. The rest can be stored at room temperature.

Procedures

1. Sample preparation
 - 1) Sample preparation from FFPE block:
 - a. Use a scalpel to trim excess paraffin off the sample block.
 - b. Cut up to 8 sections, each 5–10 µm thick. Discard the first 2–3 sections if the sample surface has been exposed to air.
 - c. Place the sections immediately in a 1.5 or 2 ml microcentrifuge tube and add 1 ml of xylene to the sample. Close the lid and vortex vigorously for 10 seconds.
 - d. Centrifuge at maximum speed for 2 minutes at room temperature.
 - e. Carefully remove the supernatant without disturbing the pellets.
 - f. Add 1 ml of ethanol (96–100%) to the pellet and mix by vortexing to extract residual xylene from the sample.
 - g. Centrifuge at maximum speed for 2 minutes at room temperature.
 - h. Carefully remove the supernatant without disturbing the pellet. Remove any remaining ethanol with a fine pipette tip.
 - i. Open the tube and incubate at room temperature or up to 37°C for 10 minutes or until all residual ethanol has evaporated.
 - 2) Sample preparation from FFPE sections on slides:
 - a. Submerge the slides in xylene I for 3 minutes, followed by xylene II for an additional 3 minutes.
 - b. Remove xylene by rinsing with 100% ethanol (1 minutes each, repeated twice).
 - c. Air dry the slides for 3-5 minutes.
 - d. Gently detach the tissue sections from the slides using a small blade, then transfer the tissue pellets into a 1.5 ml microcentrifuge tube.
2. Proteinase K solution preparation: Combine 260 µl of Proteinase K Buffer with 5 mg of Proteinase K Powder. Vortex the mixture to ensure complete dissolution. Store the solution at -20°C.
3. Resuspend the pellets in 180 µl Buffer VTL. Add 20 µl of proteinase K solution and mix by vortexing.
4. Incubate at 56°C for 1 hour or until the sample is completely lysed.
5. Incubate at 90°C for 1 hour without agitation.
6. Cool to room temperature and briefly centrifuge the tube. For RNA-free genomic DNA, add 2 µl of RNase A (100 mg/ml) and incubate for 2 minutes at room temperature.
7. Add 200 µl of Buffer VL and 200 µl ethanol to the sample. Mix thoroughly by vortexing.
8. Carefully transfer the entire lysate to a DNA Elute column and centrifuge at 8000×g (or 10000 rpm) for 1 minute. Discard the flow-through.
9. Add 500 µl of Buffer VW1 and centrifuge at 8000×g (or 10000 rpm) for 1 minute. Discard the flow-through.
10. Add 500 µl of Buffer VW2 and centrifuge at 8000×g (or 10000 rpm) for 1 minute.

- Discard the flow-through and collection tube.
11. Place the DNA Elute column in a clean 2 ml collection tube. Centrifuge at maximum speed for 3 minutes with the lid open to completely dry the membrane.
 12. Place the DNA Elute column in a clean 1.5 ml microcentrifuge tube and apply 20–100 μ l of Buffer VTE to the center of the membrane. Ensure that Buffer VTE is at room temperature.

13. Incubate at room temperature for 5 minutes and centrifuge at maximum speed (20,000×g or 14,000 rpm) for 1 minute to elute the DNA.

Note: This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Precautions: Handle with care. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.