



VitroView™ Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit

for Osteoclast Detection (50-100 Slides)

SKU: VB-3047

Introduction

Tartrate-Resistant Acid Phosphatase (TRAP) is an acid phosphatase enzyme that remains active in the presence of tartrate. The enzyme hydrolyzes a chromogenic substrate, producing a colored reaction product at sites of TRAP activity.

TRAP-positive cells typically appear as red, dark red, or purple red, depending on the kit formulation and counterstain used. Osteoclasts are generally identified as multinucleated TRAP-positive cells.

This VitroView™ Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit is designed for histochemical identification of osteoclasts and other TRAP-positive cells in tissue sections and cell preparations.

Application

- Identification of osteoclasts in bone tissue
- Assessment of osteoclast differentiation in cell cultures
- Research applications involving bone remodeling and resorption
- Histological and cytochemical investigations

Key Advantages

- Dual-Modality Imaging: Enables both bright-field and fluorescence microscopy evaluation.
- Organic Mounting Compatible: Formulated with NewFuchsin to withstand alcohol dehydration and xylene clearing.
- High Yield Capacity: Single kit processes 50 to 100 slides efficiently.
- Versatile Sample Compatibility: Validated for cell smears, adherent cells, frozen sections, and paraffin sections.
- Streamlined Working Protocol: Micro-volume master mix preparation takes less than two minutes. Clear Cellular Contrast: Delivers distinct wine-red cytoplasm against light blue nuclei.

Kit Contents

VB-3047-1	NewFuchsin Solution	0.25 ml
VB-3047-2	Sodium Nitrite Solution	0.25 ml

VB-3047-3	Naphthol AS-BI Solution	0.25 ml
VB-3047-4	TRAP Buffer	30 ml
VB-3047-5	Mayer's Hematoxylin Solution	30 ml

Storage:

Store TRAP Buffer and Mayer's Hematoxylin solution at room temperature away from light. Store other reagents at -20°C away from light. This kit is stable for at least 3 months.

Procedure**1. Sample Preparation:**

- For Cell Smears (Blood / Bone Marrow): Prepare smears using fresh samples according to routine operations. Fix in 10% Neutral Buffered Formalin (NBF) for 15–30 minutes. Wash 3 times with distilled water. Proceed to Procedure Step 2.
- For Adherent Cells / Coverslips: Discard the culture medium completely. Wash 3–4 times with PBS. Fix in 10% Neutral Buffered Formalin (NBF) for 15–30 minutes. Wash 3 times with distilled water. Proceed to Procedure Step 2.
- For Frozen Sections: Warm the frozen sections to room temperature. Fix in 10% Neutral Buffered Formalin (NBF) for 15–30 minutes. Wash 3 times with distilled water. Proceed to Procedure Step 2.
- For Paraffin Sections: Deparaffinize sections in xylene (2 × 6 min), followed by rehydration in 100% ethanol (2 min), 95% ethanol (2 × 2 min), and 70% ethanol (2 min). Rinse in distilled water for 5 min and proceed to Procedure Step 2.

2. Staining

- 1) Prepare approximately 1.0 mL of New Fuchsin TRAP staining working solution, sufficient for staining 2-5 slides, as follows:
 - In a microcentrifuge tube, combine 10 µL of NewFuchsin Solution with 10 µL of Sodium Nitrite Solution. Incubate the mixture at room temperature for 1 min.
 - Add 1.0 mL of TRAP Buffer to the tube.
 - Add 10 µL of Naphthol AS-BI Solution and mix thoroughly to prepare the working staining solution.
- 2) Use hydrophobic barrier pen to draw a water-repellent circle around tissue sections or cells on the slide.
- 3) Gently drop the working solution to cover the cells or tissue section on the glass slides and incubate at 37 °C in moisture chamber for 40–60 min.
- 4) Drop Mayer's Hematoxylin Solution onto the bone sections for 2–5 min; then wash the samples with running water for 15 min.

3. Dehydration and mounting

- 1) Dehydrate with 2 changes of 95% Ethanol and 2 changes of 100% Ethanol (2 minutes per change).
- 2) Clear with 3 changes of xylene (5 minutes per change)
- 3) Mount coverslip onto glass slide with Permount or some other suitable organic mounting medium.

4. Observation: Bright-field Microscopy can be used to examine specimens. When observing fluorescence, use a rhodamine excitation filter (500–570 nm).

Expected Results under Bright-field Microscopy

- Osteoclast cytoplasm ----wine-red
- Nuclei----- light blue

Expected Results under Fluorescence Microscopy

- Osteoclast cytoplasm ----- red

Positive Controls

- Mouse fetus whole sections (Spinal bone)
- Metaphyses or growth plates from juvenile mice or rats (3 to 6 weeks old) are the most common laboratory controls.
- Human giant cell tumor tissue sections

References

1. Nakamura A, et al (2025). Osteoclast visualization: Tartrate-resistant acid phosphatase activity staining using NewFuchsin compatible with non-aqueous mounting and tissue clearing. *Methods X*, 14: 103136
2. Luo G, et al (2025). Precision-targeting and dual silencing osteoclastogenesis and inflammatory pathways for the treatment of radiation-induced bone deterioration. *Biomaterials Advances*, 117:214369.

Note: This product is intended for research purposes only. This product is **not** intended to be used for therapeutic or diagnostic purposes in humans or animals.

Precautions: Handle with care. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.