



Data Sheet

Product Name: BrdU Positive Control FFPE Slides

SKU#: MsPs-BP-01

Product type: FFPE Sections

Species: Mouse

Slide Preparation: FVB mouse was injected with BrdU via intraperitoneal injection. Mouse intestine tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μ m in thickness and mounted each on positively charged glass slides.

Quality Control:

At least one of the tissue slides from each lot was stained with in situ BrdU IHC to ensure the quality

Package size

- 1 section per slide
- 5 slides per package

Applications

In situ BrdU positive control

Storage

2-8°C

Note: This product is intended for research purposes only.