



### **VitroView™ PBS (10×) 1L**

**SKU#: VB-6000**

#### **Description**

The packet contains a 10× concentrated formulation of phosphate buffered saline (PBS) that can be used as a wash buffer during immuno-labeling procedures.

#### **Contents**

PBS (10×), pH 7.4-----1L

#### **Storage**

Store at room temperature

#### **Procedure**

This is a 10× stock solution and should be diluted 10-fold with distilled water to make the 1× working solution.

**Note:** This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.