



VitroView™ TBS (10×), 1L

SKU#: VB-6001

Description

The packet contains a 10× concentrated formulation of Tris buffered saline (TBS) that can be used as a wash buffer during immuno-labeling procedures.

Contents

TBS (10×), pH 7.4 1L

Storage

Store at room temperature

Procedure

This is a 10× stock solution and should be diluted 10-fold with distilled water to make the 1× working solution.

Note: This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.