



VitroView™ Antibody Dilution Buffer (100ml)

SKU#: VB-6002

Description:

This antibody dilution buffer is ready-to-use solution for dilution of primary and secondary antibodies in immunoassays including immunofluorescence(IF), immunohistochemistry (IHC) and Western blot(WB).

Components

- BSA (IHC grade, IgG-Free, Protease -free)
- Cold fish skin gelatin
- Sodium azide
- 0.01M PBS, pH7.2

Storage

Store at 2-8°C.

Application:

IHC, ICC/IF, WB. The usefulness of this product in other applications has not been determined.

Note:

1. The sodium azide is an inhibitor of HRP. **Do not** use it to dilute HRP conjugated antibody.
2. The phosphate in PBS is an inhibitor of Alkaline Phosphatase (AP). **Do not** use it to dilute AP conjugated antibody.

Caution: This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.