



VitroView™ Antibody Dilution Buffer (100ml)

SKU#: VB-6002

Description:

This antibody dilution buffer is ready-to-use solution for dilution of primary and secondary antibodies in immunoassays including immunofluorescence(IF), immunohistochemistry (IHC) and Western blot(WB).

Components

- BSA (IHC grade, IgG-Free, Protease -free)
- Cold fish skin gelatin
- Sodium azide
- 0.01M PBS, pH7.2

Storage

Store at 2-8°C.

Application:

IHC, ICC/IF, WB. The usefulness of this product in other applications has not been determined.

Note:

1. The sodium azide is an inhibitor of HRP. **Do not** use it to dilute HRP conjugated antibody.
2. The phosphate in PBS is an inhibitor of Alkaline Phosphatase (AP). **Do not** use it to dilute AP conjugated antibody.

Caution: This product is intended for research purposes only. This product is **not** intended to be used for therapeutic or diagnostic purposes in humans or animals. Avoid contact with eyes, skin and clothing. Do not ingest. Wear gloves.