



Data Sheet

Product Name: Human Breast Cancer with ER+/PR+/Her2 (-) & Normal Matched Pair Tissue FFPE Sections
SKU#: HuPS-02001TA

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	HuPS-02001T	Human breast invasive ductal carcinoma and adjacent normal tissue	Female/45	2	ER(80% +), PR(70%+), E-cad(+), P120(+), CK5/6(-), CK8/18(+), Ki67(15% +), Her2 (-)
Normal Matched Tissue	HuPS-02001A	Adjacent normal tissue			

Package Size: 2 sections per slide and 2x5 slides per package
 Breast Cancer (HuPS-02001TA) ----- 5 slides
 Normal Matched Tissue (HuPS-02001A) -----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.