



Data Sheet

Product Name: Human Breast Cancer with ER (-)/PR (-)/Her2 (1+)/p53 (++) & Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-02002TA

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	TMN	IHC Data
Breast Cancer	HuPS-02002TA	Human breast invasive ductal carcinoma	Female/37	3	N/A	ER(-), PR(-), E-cad(+), P120(+), CK5/6(++), VEGF(-), Ki67(55% +), Her2 (1+).p53 (++)
Normal Matched Tissue	HuPS-02002A	Adjacent normal tissue				

Package size: 2 sections per slide, 2x5 slides per package
 Breast Cancer (HuPS-02002TA) ----- 5 slides
 Normal Matched Tissue (HuPS-02002A) -----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.