



Data Sheet

Product Name: Human Breast Cancer, ER (+)/PR(1+)/P53(+), FFPE Sections

SKU#: HuPS-02004T

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	Human breast Invasive ductal carcinoma	Female/50	1	ER(100% +), PR(100%+), P53(+), CK5/6(-),P63(-), SMA(-),Ki67(25% +),Her2 (+-),VEGF(+),E-cad(+)

Package size

- 2 sections per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.