



## Data Sheet

**Product Name:** Human Breast Cancer with ER+/PR+/Her2 (+) & Normal Matched Pair Tissue FFPE Sections

**SKU#:** HuPS-02005TA

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Breast

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	HuPS-02005TA	Human breast invasive ductal carcinoma and adjacent normal tissue	Female/50	2	ER(30% +), PR(30%+), CK-7(+),P53(-), CK5/6(-), Ki67(25% +),Her2 (++) , VEGF(+), E-cad(+)
Normal Matched Tissue	HuPS-02005A	Adjacent normal tissue			

**Package Size:** 2 sections per slide, 2x5 slides per package  
 Breast Cancer (HuPS-02005T) ----- 5 slides  
 Normal Matched Tissue (HuPS-02005A) -----5 slides

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.