



## Data Sheet

**Product Name:** Human Lung Adenocarcinoma FFPE Sections

**SKU#:** HuPS-03006T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Lung

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5  $\mu$ m in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

| Tissue      | Pathology Diagnosis | Gender/Age | IHC Data |
|-------------|---------------------|------------|----------|
| Lung Cancer | Lung adenocarcinoma | Male/53    | Ki67     |

**Package size**

- 2 sections per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.