



Data Sheet

Product Name: Human Liver Cancer FFPE Sections

SKU#: HuPS-08001T

Product type: FFPE Sections

Species: Human

Tissue Name: Liver

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

| Tissue | Pathology Diagnosis | Gender/Age | IHC Data |
|--------------|--------------------------|------------|----------|
| Liver Cancer | Hepatocellular carcinoma | Male/46 | Ki67 |

Package size

- 2 sections per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.