



## Data Sheet

**Product Name:** Human Liver Cancer & Normal Matched Pair Tissue FFPE Sections

**SKU#:** HuPS-08002TA

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Liver

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Part No	Pathology Diagnosis	Gender/Age	IHC Data
Liver Cancer	HuPS-08002T	Hepatocellular carcinoma	Male/42	Ki67
Normal Matched Tissue	HuPS-08002A	Adjacent normal tissue		

**Components:** 2x5 slides per package

Liver Cancer (HuPS-08002T) ----- 5 slides

Normal Matched Tissue (HuPS-08002A) -----5 slides

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.