



Data Sheet

Product Name: Human Liver Cancer & Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-08004TA

Product type: FFPE Sections

Species: Human

Tissue Name: Liver

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

| Tissue | Part No | Pathology Diagnosis | Gender/Age | Grade | TMN | IHC Data |
|-----------------------|--------------|---|------------|-------|-----|----------|
| Liver Cancer | HuPS-08004TA | Hepatocellular carcinoma and adjacent normal tissue | Male/56 | | | Ki67 |
| Normal Matched Tissue | HuPS-02001A | Adjacent normal tissue | | | | |

Components: 2x5 slides per package

Liver Cancer (HuPS-08004TA) ----- 5 slides

Normal Matched Tissue (HuPS-02001A) -----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.