



**Data Sheet**

Product Name: Mouse Whole Embryo (E13.5) Paraffin Sagittal Sections (1 section/slide, 5 slides/package)

SKU#: MsPs-E13.5-1A

Product type: FFPE Sections

Species: FVB Mouse

Tissue Name: Embryo (E13.5)

**Slide Preparation:** Mouse tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

**Donor Information:**

Tissue	Pathology Diagnosis	Age	IHC Data
FVB Mouse Embryo (E13.5)	Normal mouse embryo	E13.5 (Day 13.5)	Ki67

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.