



Data Sheet

Product Name: VitroView™ DNase-Treated TUNEL Positive Control Slides

SKU#: MsPs-TP-01

Product type: FFPE Sections

Species: Mouse

Slide Preparation: Mouse tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides. Mouse FFPE tissue sections were treated with DNase I. Following DNase I treatment, the slides were dehydrated and covered with paraffin for long term storage.

Quality Control: At least one of the tissue slides from each lot was stained with in situ TUNEL assay to ensure the quality.

Donor Information: mouse

Package size

- 1 section per slide
- 5 slides per package

Applications

In situ TUNEL Positive Control Slides

Storage

2-8°C

Note: This product is intended for research purposes only.