



Data Sheet

Product Name: Human Breast Cancer, ER+/PR+/Her2 (-), FFPE Sections

SKU#: HuPS-02001T1

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Grade	TMN	IHC Data
Human Breast	Human breast Invasive ductal carcinoma	Female/45	2	N/A	ER(80%+), PR(70%+), E-cad(+), P120(+), CK5/6(-),CK8/18(+), Ki67(15%+),Her2 (-)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.