



## Data Sheet

**Product Name:** Human Breast Cancer, ER(+)/PR(+)/Her2 (-) P53(+), FFPE Sections

**SKU#:** HuPS-02003T1

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Breast

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Gender/Age (year)	Tissue of Origin	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Breast Cancer	Human breast Invasive ductal carcinoma	Female/43	Human right breast	30%	I	N/A	ER(95% +), PR(+++), P53(+) CK5/6(-),P63(-), Ki67(30% +),Her2 (-),VEGF(+),P120(+)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.