



## Data Sheet

**Product Name:** Human Breast Cancer with ER(+)/PR(1+)/P53(+)& Normal Matched Pair Tissue FFPE Sections  
**SKU#:** HuPS-02004TA

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Breast

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	HuPS-02004T	Human breast Invasive ductal carcinoma	Female/50	1	ER(100%+), PR(100%+), P53(+), CK5/6(-), P63(-), SMA(-), Ki67(25%+), Her2 (+-), VEGF(+), E-cad(+)
Normal Matched Tissue	HuPS-02004A	Adjacent normal tissue			

**Package Size:** 1 section per slide, 2x5 slides per package

Breast Cancer (HuPS-02004T)-----5 slides

Normal Matched Tissue (HuPS-02004A)----- 5 slides

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.