



Data Sheet

Product Name: Human Breast Cancer with ER (+)/PR (1+)/HER2 (-)/P53 (+) & Normal Matched Pair Tissue FFPE Sections
SKU#: HuPS-02004TA

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	HuPS-02004T	Human breast Invasive ductal carcinoma	Female/50	1	ER(100% +), PR(100%+), P53(+), CK5/6(-),P63(-), SMA(-),Ki67(25% +),Her2 (-),VEGF(+),E-cad(+)
Normal Matched Tissue	HuPS-02004A	Adjacent normal tissue			

Package Size: 1 section per slide, 2×5 slides per package
 Breast Cancer (HuPS-02004T) ----- 5 slides
 Normal Matched Tissue (HuPS-02004A) -----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.