



Data Sheet

Product Name: Human Breast Cancer, ER+/PR+/Her2(+), FFPE Sections

SKU#: HuPS-02005T1

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	Human breast Invasive ductal carcinoma	Female/50	2	ER(30% +), PR(30%+), CK-7(+),P53(-), CK5/6(-), Ki67(25% +),Her2 (++) , VEGF(+), E-cad(+)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.