



Data Sheet

Product Name: Human Breast Cancer with ER+/PR+/Her2(+)& Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-02005TA

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Breast Cancer	HuPS-02005TA	Human breast invasive ductal carcinoma and adjacent normal tissue	Female/50	2	ER(30% +), PR(30%+), CK-7(+),P53(-), CK5/6(-), Ki67(25% +),Her2 (++) , VEGF(+), E-cad(+)
Normal Matched Tissue	HuPS-02005A	Adjacent normal tissue			

Package Size: 1 section per slide, 2x5 slides per package

Breast Cancer (HuPS-02005T)-----5 slides

Normal Matched Tissue (HuPS-02005A)----- 5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.