



## Data Sheet

**Product Name:** Human Lung Adenocarcinoma FFPE Sections

**SKU#:** HuPS-03001T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Lung

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Organ& Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Biomarker
Lung Cancer	Adenocarcinoma	Female/73			2	Ki67(+) 60%, CK7 (+), EGFR(+),TTF-1 (+)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH, ISH, RNAscope ISH and DNA/RNA isolation for NGS

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.