



Data Sheet

Product Name: Human Lung Adenocarcinoma FFPE Sections

SKU#: HuPS-03003T1

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μ m in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Smoking	Grade	TMN/Stage	Biomarker
Lung	Human Lung Adenocarcinoma	Male/51		2		Ki67(+, 15%), CK (+), TTF-1 (+), P40 (-), CK5/6(-)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.