



Data Sheet

Product Name: Human Lung Adenocarcinoma FFPE Sections

SKU#: HuPS-03004T

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μ m in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ& Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Biomarker
Lung Cancer	Adenocarcinoma	Male/57			2	CK (+), CK7(+),CK20 (-), TTF-1 (+), Ki67(+)

Package size

- 1 sections per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH, ISH and NGS, etc.

Storage

2-8°C

Note: This product is intended for research purposes only.