



Data Sheet

Product Name: Human Lung Squamous Carcinoma FFPE Sections

SKU#: HuPS-03005T

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ& Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Biomarker
Lung Cancer	Lung squamous carcinoma	Male/57			2	CK (+), P63(+), CK5/7(+), TTF-1 (-), Napsin A (-), CD56 (-), Ki67(+)

Package size

- 1 sections per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.