



## Data Sheet

**Product Name:** Human Lung Squamous Carcinoma FFPE Sections

**SKU#:** HuPS-03005T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Lung

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Organ&Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Biomarker
Lung Cancer	Lung squamous carcinoma	Male/57	NA	NA	2	CK (+),P63(+), CK5/7(+), TTF-1 (-), Napsin A (-), CD56 (-), Ki67(+)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.