



Data Sheet

Product Name: Human Lung Squamous Carcinoma & Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-03005TA

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	IHC Data
Lung Cancer	HuPS-03005T	lung squamous carcinoma and adjacent non-cancer lung tissue	Male/57	Ki67
Normal Matched Tissue	HuPS-03005A	Adjacent normal tissue		

Package Size: 2 sections per slide, 2x5 slides per package
 Lung Cancer (HuPS-03005T) ----- 5 slides
 Normal Matched Tissue (HuPS-03005A) -----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.