



### Data Sheet

**Product Name:** Human Lung Squamous Carcinoma & Adjacent Matched Pair Non-Tumor (Chronic Inflammation) Tissue FFPE Sections

**SKU#:** HuPS-03005TA

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Lung

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Part No	Pathology Diagnosis	Gender/Age	IHC Data
Lung Cancer	HuPS-03005T	lung squamous carcinoma and adjacent non-cancer lung tissue	Male/57	Ki67
Adjacent Matched non-tumor Tissue	HuPS-03005A	Adjacent non-tumor (chronic inflammation) Tissue		

**Package Size:** 2 sections per slide, 2x5 slides per package

Lung Cancer (HuPS-03005T) ----- 5 slides

Non-tumor Matched Non-Tumor Tissue (HuPS-03005A) -----5 slides

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.