



Data Sheet

Product Name: Human Colon Mucinous Adenocarcinoma & Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-07002TA

Product type: FFPE Sections

Species: Human

Tissue Name: Colon

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	SKU No	Pathology Diagnosis	Gender/Age	IHC Data
Colon Cancer	HuPS-07002T	Colon mucinous adenocarcinoma	Male/57	Ki67
Normal Matched Tissue	HuPS-07002A	Adjacent normal tissue		

Package Size: 1 section per slide, 2x5 slides per package

Colon Cancer (HuPS-07002T)-----5 slides
 Normal Matched Tissue (HuPS-07002A)-----5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.