



## Data Sheet

**Product Name:** Human Colon Adenocarcinoma & Normal Matched Pair Tissue FFPE Sections

**SKU#:** HuPS-07005TA

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Colon

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Part No	Pathology Diagnosis	Gender/Age	IHC Data
Colon Cancer	HuPS-07005T	Colon adenocarcinoma and adjacent normal tissue	Female/51	Ki67
Normal Matched Tissue	HuPS-07005A	Adjacent normal tissue		

**Package Size:** 1 section per slide, 2x5 slides per package

Colon Cancer (HuPS-07005TA)-----5 slides

Normal Matched Tissue (HuPS-07005A)-----5 slides

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.