



## Data Sheet

**Product Name:** Human Liver Cancer FFPE Sections

**SKU#:** HuPS-08002T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Liver

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5  $\mu\text{m}$  in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Tissue	Pathology Diagnosis	Gender/Age	IHC Data
Liver Cancer	Hepatocellular carcinoma	Male/42	Ki67

**Package size**

- 2 sections per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.