



Data Sheet

Product Name: Human Liver Cancer & Normal Matched Pair Tissue FFPE Sections

SKU#: HuPS-08003TA

Product type: FFPE Sections

Species: Human

Tissue Name: Liver

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Part No	Pathology Diagnosis	Gender/Age	Grade	IHC Data
Liver Cancer	HuPS-08003T	Hepatocellular carcinoma	Female/65	2	Ki67
Normal Matched Tissue	HuPS-02001A	Adjacent normal tissue			

Package Size: 1 section per slide, 2x5 slides per package

Liver Cancer (HuPS-08003T)-----5 slides

Normal Matched Tissue (HuPS-08003A)----- 5 slides

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.