



Data Sheet

Product Name: Mouse Whole Embryo (E18.5) Paraffin Sagittal Sections (1 section/slide, 5 slides/package)

SKU#: MsPs-E18.5-3A

Product type: FFPE Sections

Species: FVB Mouse

Tissue Name: Embryo (E18.5)

Slide Preparation: Mouse tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Age	IHC Data
FVB Mouse Embryo (E18.5)	Normal mouse embryo	E18.5 (Day 18.5)	Ki67

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.