



**Data Sheet**

Product Name: Human Skin Tissue (Normal) FFPE Sections

SKU#: HuPS-01001A1

Product type: FFPE Sections

Species: Human

Tissue Name: Human Abdomen Skin

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

**Donor Information:**

Tissue	Pathology Diagnosis	Gender/Age	IHC Data
Human Skin	Normal skin tissue	Female/39	CK14

**Package size**

- 1 Section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.