



**Data Sheet**

**Product Name:** Human Skin Tissue (Normal) FFPE Block

**SKU#:** HuPS-01001A-B

**Product type:** FFPE Block

**Species:** Human

**Tissue Name:** Human Abdomen Skin

**FFPE Block Preparation:** Tissue is prepared by a certified histotechnologist under supervision of a medical pathologist. FFPE block preparation: Tissues are excised, immediately fixed in 10% neutral formalin for 24 hours, processed in a automatic tissue processor and embedded in paraffin. The target FFPE specimen size is 0.5 x 1 x 1 cm, but it can vary significantly depending on the organ.

**Quality Control:**

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

**Donor Information:**

Tissue	Pathology Diagnosis	Gender/Age	IHC Data
Human Skin	Normal skin tissue	Female/39	CK14

**Package size**

1 FFPE block per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.