



Data Sheet

Product Name: Human Umbilical Cord Tissue (Normal) FFPE Block

SKU#: HuPS-10001-B

Product type: FFPE Block

Species: Human

Tissue Name: Human umbilical cord

FFPE Block Preparation: Tissue is prepared by a certified histotechnologist under supervision of a medical pathologist. FFPE block preparation: Tissues are excised, immediately fixed in 10% neutral formalin for 24 hours, processed in a automatic tissue processor and embedded in paraffin. The target FFPE specimen size is 0.5 x 1 x 1 cm, but it can vary significantly depending on the organ.

Quality Control:

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Gestation (Weeks)	IHC test
Human umbilical cord	Normal umbilical cord tissue	Female/22	40	yes

Package size

1 FFPE block per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.