



Data Sheet

Product Name: Human Umbilical Cord Tissue (Normal) FFPE Sections

SKU#: HuPS-11001

Product type: FFPE Sections

Species: Human

Tissue Name: Umbilical cord

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	Race	Weeks Gestation	IHC Data
Human umbilical cord	Normal umbilical cord tissue	Male/22YO	Caucasian	40	Yes

Package size

- 1 Section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.