



**Data Sheet**

Product Name: Human Placenta Tissue (Normal) FFPE Sections

SKU#: HuPS-11001

Product type: FFPE Sections

Species: Human

Tissue Name: Placenta

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

**Donor Information:**

Tissue	Pathology Diagnosis	Gender/Age	Race	Week Gestation	IHC Data
Human Placenta	Normal placenta tissue	Male/21YO	Caucasian	39.3	

**Package size**

- 1 Section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.