



Data Sheet

Product Name: Human Lung Adenocarcinoma FFPE Sections

SKU#: HuPS-03006T

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Tested Biomarker
Lung Cancer	Lung adenocarcinoma	Male/53			2	CK7 (+), CK20 (-), CK5/6(-), P63 (-), TTF-1 (+), Napsin A(-), Ki 67 (+, 35%)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH, ISH and NGS

Storage

2-8°C

Note: This product is intended for research purposes only.