



Data Sheet

Product Name: Human Breast Cancer, Triple Positive (ER+/PR+/Her2+), FFPE Sections

SKU#: HuPS-02014T

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Therapies (Chemotherapies ,radiation and Hormone) before Collection of Sample	Gender/Age (year)	Tumor area	Tumor Grade	TMN Stage	ER,PR & HER2 IHC Data
Breast	Multifocal infiltrating duct carcinoma	No	Female/50	50%	II	T4N2 M0	ER(95% +++), PR(50%+++), Her2 (Score 2+)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.