



**Data Sheet**

Product Name: Human Breast Cancer, Triple Positive (ER+/PR+/Her2+), FFPE Sections

SKU#: HuPS-02015T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Breast

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Therapies (Chemotherapies ,radiation and Hormone) before Collection of Sample	Gender/Age (year)	Tumor area	Tumor Grade	TMN Stage	ER,PR & HER2 IHC Data
Breast	Infiltrating duct carcinoma	No	Female/68	45%	III	T3pN 3aM0	ER(80% ++), PR(70%+++), Her2 (Score 2)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.