



Data Sheet

Product Name: Human Breast Cancer, Triple Negative (ER-/PR-/Her2-), FFPE Sections

SKU#: HuPS-02022T

Product type: FFPE Sections

Species: Human

Tissue Name: Breast

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Therapies (Chemotherapies ,radiation and Hormone) before Collection of Sample	Gender/Age (year)	Tumor area	Tumor Grade	TMN Stage	ER,PR & HER2 IHC Data
Breast	Infiltrating duct carcinoma	No	Female/73	86%	II	T4pN 1AM0	ER(-), PR(-), Her2 (-)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.