



## Data Sheet

**Product Name:** Human Small Cell Lung Cancer FFPE Sections

**SKU#:** HuPS-03032T2

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Lung

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5  $\mu$ m in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Organ /Tissue	Pathology Diagnosis	Tumor Area	Gender/ Age	Smoking	TMN/Stag e	Grade	Biomarker
Lung Cancer	Small cell lung cancer	90%	Male/65	N/A	T2M0N0 /IIB	N/A	N/A

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.