



Data Sheet

Product Name: Human Small Cell Lung Cancer FFPE Sections

SKU#: HuPS-03037T

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ /Tissue	Pathology Diagnosis	Tumor Area	Gender/ Age	Smoking	TMN/Stag e	Grade	Biomarker
Lung Cancer	Small cell lung cancer	20%	Female/ 48	N/A	T1bM0N3 /IIIB	N/A	N/A

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.