



## Data Sheet

**Product Name:** Human Pancreas Cancer FFPE Sections

**SKU#:** HuPS-12006T2

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Pancreas Cancer

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5  $\mu\text{m}$  in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Gender/Age (year)	Tissue of Origin	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Pancreas Cancer	Adenocarcinoma, Moderately Differentiated	Male/52	Human Pancreas	35%	II	N/A	

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.