



Data Sheet

Product Name: Human Pancreas Cancer FFPE Sections

SKU#: HuPS-12016T

Product type: FFPE Sections

Species: Human

Tissue Name: Pancreas Cancer

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Gender/Age (year)	Tissue of Origin	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Pancreas Cancer	Adenocarcinoma and Neuroendocrine Tumor, Well Differentiated	Male/67	Human Pancreas of ampula of vater	85%	I	T3N0M0/Stage III	

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.