



Data Sheet

Product Name: Human Gastric FFPE Sections

SKU#: HuPS-13004T

Product type: FFPE Sections

Species: Human

Tissue Name: Gastric

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

| Tissue | Pathology Diagnosis | Gender/Age (year) | % Tumor area | Tumor Grade | TMN Stage | IHC Data |
|---------------|--|-------------------|--------------|-------------|-----------|----------|
| Human stomach | Residual moderately differentiated adenocarcinoma of the stomach | Male/57 | 20% | II | T3N1MO | |

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.