



## Data Sheet

**Product Name:** Human Ovary Cancer FFPE Sections

**SKU#:** HuPS-14006T1

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Ovary

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

| Organ& Tissue | Pathology Diagnosis                  | Gender/Age | % Tumor Area | Grade | TMN Stage | Biomarkers |
|---------------|--------------------------------------|------------|--------------|-------|-----------|------------|
| Human Ovary   | endometrioid adenocarcinoma of ovary | Female/46  | 90%          | III   | T1bN0M0   |            |

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH, ISH, RNAscope ISH and DNA/RNA isolation for NGS

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.